



DREAM-C Operation Guide

Diagnostic Receptive Expressive Assessment of Mandarin – Comprehensive

Norm-referenced language screening for Mandarin-speaking children aged 2.5-8 years







• • • •





Log in

Run the software

Enter username and password



Make sure your device is turned on and connected to internet



Click on the Bethel DREAM software icon

the device needs to stay connected to the internet throughout the test

DETTILE	
Username	
Password	
Remember Password	Forgot your Password?
Log	in
No account	? Sign up

••••



Create an

Assessment

PART 02



Purchase an Assessment License



Click on the "Store" Icon and select "License".





Purchase Test Protocols / Reports





Redeem Discount Code





Create an Assessment--Select DREAM-C





Enter Tool

Click on Enter Tool

Your remaining number of reports will be displayed above the button. When report balance is insufficient, a prompt will remind you to repurchase.



Create an Assessment-- Create a New Test

0:43 3月21日 退出工具	周五			중 68% ■	Create a New Test
	Q.搜索			搜索	Click the "New Test" button at the bottom of the
	未完成	已完成			page to create a new test
	测试编	号:C10025 继续测	进行中		
			<u> 文に 7 書 </u>		
ream-C-Probation	n(com.portal.dream.tools)		利廷测试		



Create an Assessment--Enter Child's Age

・TestFlight 11:26 返回	2月27日周四	२ 73% 👍
选择	译已添加儿童 X 搜 搜 搜	Ŕ
	孩子的年龄 年: 0岁 ∨ 月: 0月 ∨	
	取消 提交	
ceam-S-Probation(com.	sportal.dream.tools)	

Enter Client's Age

Select child's age and click the "Submit" button



Create an Assessment--Prerequisites



Prerequisites

When creating a test, you need to confirm with the parents whether the child meets the basic requirements for the test. A test can only be created if all the criteria are met.

Click the "Submit" button to return to the test.

If a test is not completed within 10 days, it will automatically expire.



Assessment Procedure

PART 03



Assessment Procedure--Start a Test

0:44 3月21日月	周五			奈 68% ■
退出工具				
	Q、搜索		搜索	
	未完成 已完成			
	测试编 <mark>+</mark> :C10063	新建测试 开始测试	进行中 继续测试	
		新建测试		
ceam-C-Probation	(com.portal.dream.tools)	が「文王 7次」にい		

Start a Test

After filling in the child's information, return to the test homepage, find the target child, and select "Start Test" to enter the test interface.

Please note:

The Bethel DREAM App does not save child's identifiable information such as name, date of birth, or sex. For records keeping, **please keep details of the child information and test report number separately.**



Assessment Procedure--Confirm Child's Info



Confirm Child's Info

After selecting "Start Test," you will be directed to the child information confirmation page. Once confirmed, click "Next" to proceed with the test.

Please make sure to confirm the child's age. Once verified, click "Next" to start the test. Once the test completed, the child's information cannot be modified.



Assessment Procedure--Operation Buttons



Operation Buttons

To prevent children from being disturbed while selecting answers, the operation buttons in the assessment are hidden.

- The top left corner is for **Save and Exit**.
- The bottom left corner is for **Returning to the Previous Question**.
- The top right corner is for **Repeating the Audio**.
- The bottom right corner is for **Next Question**.

Note that every item has audio prompts. To adjust the **volume** of the prompts, use the **volume buttons** on iPad.

When saving and exiting, please ensure you have a stable internet connection to prevent data loss. All tests can be saved and resumed within 10 days.



Assessment Procedure- Complete the Test



Save & Exit

After completing the test, you can click "**View Report**" to check the report or click "**Return**" to go back to the children's list interface.



Assessment Report

PART 04

0



Assessment Report--View Report

After completing the test, click "view report" to check your Assessment report





Saving, Printing and Forwarding Reports



C100 PDF文	83 稿 · 1.1 MB			
AirDrop	Mail	备忘录	() 微信	G
Сору			ß	
Markup			\otimes	
Print			ē	
存储到"文件"				
Mobile I	Print - 打印			
在QQ邮箱中打开			\square	
编辑操作	乍			

Completed test reports will appear in the **Reports** section of the Bethel DREAM App.

- To save, print, email or send reports via other apps,
- 1. Click on the checkbox beside the report
- 2. Click on the Send button
- 3. A menu of options will appear which will allow you to email, print or forward the selected test reports, using your choice of apps.

Note that

- The Bethel DREAM App Stores test reports for 30 days – please ensure that all reports are saved/forwarded off within this time period
- The Bethel DREAM App does not save child's identifiable information such as name, date of birth, or sex. For records keeping, please keep details of the child information and test report number separately.



Please contact us via dreamapp@ourbethel.com if you have further questions

